

Vacuum Oven

진공건조기



Features

- The grounded circuit and NFB prevent a short circuit and over current of electricity and fuse change is not required.
- Equipped overheat protection sensor.
- Accurate temperature controls by digital PID controller.
- High temperature ceramic wool optimize insulation and heat loss.
- Tempered safety glass windows for excellent insulation and easy observation.
- Set time and stable temperature function.

Specifications


Made in Korea

Spec. / Model		VO-27	VO-64
Camber volume		27ℓ	64ℓ
Temp	Range	Max 200 °C	
	Controller	Digital PID Controller	
Size	Internal (WxDxH, mm)	300x300x300	400x400x400
	External (WxDxH, mm)	515x460x670	615x560x770
Vacuum range		1x10 ⁻³ ~76cmHg	
Safety	Electrical	Over temp limiter, ELB (Earth Leakage circuite breaker)	
Timer		Max 99hr 59min	
Material	Internal	Stainless steel #304	
	External	Steel, Powder coating	
	Shelves	AL 5052 2.0t	
	Heater	Sheath Heater	
	Insulation	Ceramic Wool 25mm	
	Door gasket	Molded heat resistance silicon rubber	
	Windows	Tempered Safety Glass, 12t	Tempered Safety Glass, 15t
Electric Requirements		220VAC. 50/60Hz	
Power Consumption		1.2KW	1.8KW

Vacuum Oven

진공건조기

Distinction



Structure

