

Digital Colony Counter

코로나카운터



Features

- By using of button and a marking pen, 2 counter mode are possible.
- This unit main power is 220V.
- Easy to read 4 digits LED digital display.
- Buzzer sounds on every counting.
- A marking pen is usable for user's convenience.
- Wider view of 125 mm wide lens.
- Bacteria can be counted by picking with stainless needle point pen.

Specifications

Made in Korea


Spec. / Model	HYC-560
Display type	4 Digital LED Display
Petri dish holder	Accepts 100mm dish (removable) 150mm (fixed)
Zero rest button	Yes
Illumination	Fluorescent lamp
Dimensions (WxDxH,mm)	310x290x100
Mode of Counting	Manual push button counting
	Needle-point probing
	Marking pen probing
Electric Requirements	220VAC. 50/60Hz

Standard accessories	
Magnifier 5" 1ea	Stainless Needle point pen 1ea
Power cord 1ea	Marking pen probe 1ea

Digital Colony Counter

코로니카운터

Distinction



Structure

